


**Product Name: Recombinant Human IL10RA Protein**  
**Catalog #: PEH3177**



## Summary

|                              |                                                                                                                                                                                                                                       |
|------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Name</b>                  | I10R1_HUMAN                                                                                                                                                                                                                           |
| <b>Purity</b>                | Greater than 85% as determined by reducing SDS-PAGE                                                                                                                                                                                   |
| <b>Endotoxin level</b>       | Please contact with the lab for this information                                                                                                                                                                                      |
| <b>Tag</b>                   | N-SUMO & C-His                                                                                                                                                                                                                        |
| <b>Accession #</b>           | Q13651                                                                                                                                                                                                                                |
| <b>Host</b>                  | E.coli                                                                                                                                                                                                                                |
| <b>Species</b>               | Human                                                                                                                                                                                                                                 |
| <b>Molecular Mass</b>        | 37 kDa                                                                                                                                                                                                                                |
| <b>Buffer</b>                | Phosphate buffered saline                                                                                                                                                                                                             |
| <b>Shipping</b>              | The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.                                                                                                                    |
| <b>Stability&amp;Storage</b> | Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.                        |
| <b>Sequence</b>              | SVWFEEFFH HILHWTPIPN QSESTCYEVA LLRYGIESWN SISNCSQTLS<br>YDLTAVTLDL YHSNGYRARV RAVDGSRHSN WVTNTNTRFSV DEVTLTVGSV<br>NLEIHNGFIL GKIQLPRPKM APANDTYESI FSHFREYEIA IRKVPGNFTF<br>THKKVKHENF SLLTSGEVGE FCVQVKPSVA SRSNKG MWSK EECISLTRQY |

## image



## Background

|                          |                                                                                                                                                                                                     |
|--------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Alternative Names</b> | Interleukin-10 receptor subunit alpha (IL-10 receptor subunit alpha) (IL-10R subunit alpha) (IL-10RA) (CDw210a) (Interleukin-10 receptor subunit 1) (IL-10R subunit 1) (IL-10R1) (CD antigen CD210) |
|--------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

## Note

**Product Name: Recombinant Human IL10RA Protein**  
**Catalog #: PEH3177**



---

For Research Use Only , Not for Diagnostic Use