


Summary

Name	MOBP_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	Q13875
Host	E.coli
Species	Human
Molecular Mass	63 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MSQKPAKEGP RLSKNQKYSE HFSIHCCPPF TFLNSKKEIV DRKYSICKSG CFYQKKEEDW ICCACQKTRT SRRAKSPQRP KQQAAPPAV VRAPAKPRSP PRSERQPRSP PRSERQPRSP PRSERQPRSP PRSERQPRPR PEVRPPPAKQ RPPQKSKQQP RSSPLRGP GA SRGGSPVKAS RFW

image



Background

Alternative Names Myelin-associated oligodendrocyte basic protein

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human MOB Protein
Catalog #: PEH3181

