


Product Name: Recombinant Human WRN Protein
Catalog #: PEH3189



Summary

Name	WRN_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q14191
Host	E.coli
Species	Human
Molecular Mass	60 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QQRKCPEWMN VQNKRCAVEE RKACVRKSVF EDDLPLEFT GSIVSYDAS DCSFLSEDIS MSLSDGDVVG FDMEWPPLYN RGKLGKVALI QLCVSESKCY LFHVSSMSVF PQGLKMLEN KAVKKAGVGI EGDQWKLLRD FDIKKNFVE LTDVANKKLK CTETWSLNSL VKHLLGKQLL KDKSIRCSNW SKFPLTEDQK LYAATDAYAG FIIYRNLEIL DDTVQRFAIN KEEIILLSDM NKQLTSISEE VMDLAKHLPH AFSKLENPRR VSILLKDISE NLYSLRRMII GSTNIETELR PSNNLNLLSF EDSTTGGVQQ KQIREHEVLI HVEDETWDPT LDHLAKHDGE DVLGNKVERK EDGFEDGVED NKLKENMERA CLMSLDITEH ELQILEQQSQ EEYLSDIAYK STEHLSPNDN ENDTSYVIES DEDLEMEMLK HLSPNDNEND TSYVIESDED LEMEMLSLE NLNSGTVEPT HSKCLKMERN LGLPTKEEEE

image



Background

Product Name: Recombinant Human WRN Protein
Catalog #: PEH3189



Alternative Names

Werner syndrome ATP-dependent helicase (EC 3.6.4.12) (DNA helicase, RecQ-like type 3) (RecQ3) (Exonuclease WRN) (EC 3.1.-.-) (RecQ protein-like 2)

Note

For Research Use Only , Not for Diagnostic Use