


Product Name: Recombinant Human DPYSL3 Protein
Catalog #: PEH3190



Summary

Name	DPYL3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q14195
Host	E.coli
Species	Human
Molecular Mass	63 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	RITSDRLLIK GGRIVNDDQS FYADIYMEDG LIKQIGDNLI VPGGVKTIEA NGKMVIPGGI DVHTHFQMPY KGMTTVDDFF QGTKAALAGG TTMIIDHVVP EPESLSTEAY EKWREWADGK SCCDYALHVD ITHWNDSVKQ EVQNLKDKG VNSFMVYMAY KDLYQVSNT ELYEFTCLGE LGAIQVHAE NGDIIAQEQT RMLEMGITGP EGHVLSRPEE LEAEAVFRAI TIASQTNCPL YVTKVMSKSA ADLISQARKK GNVVFGEPI T ASLGIDGTHY WSKNWAKAAA FVTSPPLSPD PTTPDYINSL LASGDLQLSG SAHCTFSTAQ KAIGKDNFTA IPEGTNGVEE RMSVIWDKAV ATGKMDENQF VAVTSTNAAK IFNLYPRKGR ISVGSDSDLV IWDPAVKIV SAKNHQSAAE YNIFEGMELR GAPLVVICQG KIMLEDGNLH VTQGAGRFIP CSPFSDYVYK RIKARRKMAD LHAVPRGMYD GPVFDLTTTP KGGTPAGSAR GSPTRPNPPV RNLHQSGFSL SGTQVDEGVR SASKRIVAPP

image



Background

Product Name: Recombinant Human DPYSL3 Protein
Catalog #: PEH3190



Alternative Names

Dihydropyrimidinase-related protein 3 (DRP-3) (Collapsin response mediator protein 4) (CRMP-4) (Unc-33-like phosphoprotein 1) (ULIP-1)

Note

For Research Use Only , Not for Diagnostic Use