


Product Name: Recombinant Human FAM50A Protein
Catalog #: PEH3191



Summary

Name	FA50A_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q14320
Host	E.coli
Species	Human
Molecular Mass	27 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MKQRIAEENI MKSNIDKKFS AHYDAVEAEL KSSTVGLVTL NDMKAKQEAL VKEREKQLAK KEQSKELQMK LEKLREKERK KEAKRKISSL SFTLEEEEEEG GEEEEEAAMY EEEMEREEIT TKKRKLGNP DVDTSFLPDR DREEENRLR EELRQWEAK QEKIKSEEIE ITFSYWDGSG HRRTVKMRKG NTMQQFLQKA LEILRKDFSE

image



Background

Alternative Names Protein FAM50A (Protein HXC-26) (Protein XAP-5)

Note

Product Name: Recombinant Human FAM50A Protein
Catalog #: PEH3191



For Research Use Only , Not for Diagnostic Use