


Product Name: Recombinant Human NTRK3 Protein
Catalog #: PEH3202



Summary

Name	NTRK3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	Q16288
Host	E.coli
Species	Human
Molecular Mass	55 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GNLFPLLEGQ DSGNSNGNAS INITDISRNI TSIHIENWRS LHTLNAVDME LYTGLQKLT I KNSGLRSIQP RAFAKNPHLR YINLSSNRLT TLSWQLFQTL SLRELQLEQN FFNCSCDIRW MQLWQEQGEA KLNSQNLICI NADGSQPLPF RMNISQCDLP EISVSHVNLT VREGDNAVIT CNGSGSPLPD VDWIVTGLQS INTHQTNLNW TNVHAINLTL VNVTSEDNGF TLTCIAENVV GMSNASVALT VYYPPRVVSL EEPELRLEHC IEFVVRGNPP PTLHWLHNGQ PLRESKIIHV EYYQEGEISE GCLLFNKPTH YNNGNYTLIA KNPLGTANQT INGHFLKEPF PESTDNFILF DEVSPTPPIT

image



Background

Alternative Names NT-3 growth factor receptor (EC 2.7.10.1) (GP145-TrkC) (Trk-C)

Product Name: Recombinant Human NTRK3 Protein
Catalog #: PEH3202



(Neurotrophic tyrosine kinase receptor type 3) (TrkC tyrosine kinase)

Note

For Research Use Only , Not for Diagnostic Use