


Summary

Name	MAPK3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q16644
Host	E.coli
Species	Human
Molecular Mass	44 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MDGETAEEQG GPVPPVAPG GPGLGGAPGG RREPKKYAVT DDYQLSKQVL GLGVNGKVLE CFHRTGQKC ALKLLYDSPK ARQVDHHWQ ASGGPHIVCI LDVYENMHHG KRCLLIIMEC MEGGELFSRI QERGDQAFTE REAAEIMRDI GTAIQFLHSH NIAHRDVKPE NLLYTSKEKD AVLKLTDFGF AKETTQNALQ TPCYTPYYVA PEVLGPEKYD KSCDMWSLGV IMYILLCGFP PFYSNTGQAI SPGMKRRIRL GQYGFPNPEW SEVSEDAKQL IRLLLKTDPT ERLTITQFMN HPWINQSMVV PQTPLHTARV LQEDKDHWE VKEEMTSALA TMRVDYDQVK IKDLKTSNNR LLNKRRKKQA GSSASQGCN NQ

image



Background

Alternative Names MAP kinase-activated protein kinase 3 (MAPK-activated protein kinase 3)



Product Name: Recombinant Human MAPKAPK3 Protein
Catalog #: PEH3208

(MAPKAP kinase 3) (MAPKAP-K3) (MAPKAPK-3) (MK-3) (EC 2.7.11.1)
(Chromosome 3p kinase) (3pK)

Note

For Research Use Only , Not for Diagnostic Use