


Product Name: Recombinant Human CSPG4 Protein
Catalog #: PEH3214



Summary

Name	CSPG4_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q6UVK1
Host	E.coli
Species	Human
Molecular Mass	10 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	TLTTAHLEAT LEEAGPSPPT FHYEVVQAPR KGNLQLQGTR LSDGQGFTQD DIQAGRVTYG ATARASEAVE

image



Background

Alternative Names	Chondroitin sulfate proteoglycan 4 (Chondroitin sulfate proteoglycan NG2) (Melanoma chondroitin sulfate proteoglycan) (Melanoma-associated chondroitin sulfate proteoglycan)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use