


Product Name: Recombinant Human SPAG16 Protein
Catalog #: PEH3224



Summary

Name	SPG16_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8N0X2
Host	E.coli
Species	Human
Molecular Mass	26 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	AVRVLEEALG MGLTAAGDAR DTADAVAAEG AYYLEQVTIT EASEDDYEYE EIPDDNFSIP EGEDLAKAI QMAQEATDT EILERKTVLP SKHAVPEVIE DFLCNFLIKM GMTRTLCFQ SEWYELIQKG VTELRTVGNV PDVYTQIMLL ENENKNLKGD LKHYKQADK AREDLLKIQK ERDFHRMHHK RIVQEKNKLI NDLKGLKLHY

image



Background

Alternative Names Sperm-associated antigen 16 protein (Pf20 protein homolog)

Note

Product Name: Recombinant Human SPAG16 Protein
Catalog #: PEH3224



For Research Use Only , Not for Diagnostic Use