


**Product Name: Recombinant Human PDIK1L Protein**  
**Catalog #: PEH3225**



## Summary

<b>Name</b>	PDK1L_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-His
<b>Accession #</b>	Q8N165
<b>Host</b>	E.coli
<b>Species</b>	Human
<b>Molecular Mass</b>	39 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	IREVGRGSYG VVYEAVIRKT SARVAVKKIR CHAPENVELA LREFWALSSI KSQHNPVIHL EECILQKDGM VQKMSHGSNS SLYLQLVETS LKGEIAFDPR SAYYLWVMD FCDGGDMNEY LLSRKPNRKT NTSFMLQLSS ALAFLHKNQI IHRDLKPDNI LISQTRLDS DLEPTLKVAD FGLSKVCSAS GQNPEEPVSV NKCFLSTACG TDFYMAPEVW EGHYAKADI FALGIIWAM LERITFIDTE TKKELLGSYV KQGTEIVPVG EALLENPKME LLIPVKKKSM NGRMKQLIKE MLAANPQDRP DAFELELRLV QIAFKDSSWE

## image



## Background

**Alternative Names** Serine/threonine-protein kinase PDIK1L (EC 2.7.11.1) (PDLIM1-interacting kinase 1-like)

**Product Name: Recombinant Human PDIK1L Protein**  
**Catalog #: PEH3225**

---



**Note**

For Research Use Only , Not for Diagnostic Use