



Product Name: Recombinant Human DPP10 Protein
Catalog #: PEH3227



Summary

Name	DPP10_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	Q8N608
Host	E.coli
Species	Human
Molecular Mass	48 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LFCEGPRVPV VSLHSTDNPA KYFILESNSM LKEAILKKKI GKPEIKILHI DDYELPLQLS LPKDFMDRNQ YALLIMDEE PGGQLVTDKF HIDWDSVLID MDNVIVARFD GRGSGFQGLK ILQEIHRRLG SVEVKDQITA VKFLLKLPYI DSKRLSIFGK GYGGYIASMI LKSDEKLFKC GSVVAPITDL KLYASAFSER YLGMPskees TYQAASVLHN VHGLKEENIL IIHGTAATKV HFQHSaelik HLIKAGVNYT MQVYPDEGHN VSEKSKYHLY STILKFFSDC LKEEISVLPQ

image



Background

Alternative Names	Inactive dipeptidyl peptidase 10 (Dipeptidyl peptidase IV-related protein 3) (DPRP-3) (Dipeptidyl peptidase X) (DPP X) (Dipeptidyl peptidase-like protein 2) (DPL2)
--------------------------	---

Product Name: Recombinant Human DPP10 Protein
Catalog #: PEH3227



Note

For Research Use Only , Not for Diagnostic Use