


Product Name: Recombinant Human PNMA1 Protein
Catalog #: PEH3232



Summary

Name	PNMA1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8ND90
Host	E.coli
Species	Human
Molecular Mass	21 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GPGEETFDPW LEHTNEVLEE WQVSDVEKRR RLMESLRGPA ADVIRILKSN NPAITTAECL KALEQVFGSV ESSRDAQIKF LNTYQNPGEK LSAYVIRLEP LLQKVVEKGA IDKDNVNQAR LEQVIAGANH SGAIRRQLWL TGAGEGPAPN LFQLLVQIRE EEAKEEEEEEA

image



Background

Alternative Names Paraneoplastic antigen Ma1 (37 kDa neuronal protein) (Neuron- and testis-specific protein 1)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human PNMA1 Protein
Catalog #: PEH3232

