



Product Name: Recombinant Human DNER Protein
Catalog #: PEH3233



Summary

Name	DNER_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	Q8NFT8
Host	E.coli
Species	Human
Molecular Mass	38 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	RQLQVPATQ EPDKILPRSQ ATVTLPWQP KTGQKVVEMK WDQVEVIPDI ACGNASSNSS AGGRLVSFEV PQNTSVKIRQ DATASLILLW KVTATGFQQC SLIDGRSVTP LQASGGLVLL EEMLALGNNH FIGFVNDSTV KSIVALRLTL VVKVSTCVPG ESHANDLECS GKGKCTTKPS EATFSCTCEE QYVGTFCEEY DACQRKPCQN NASCIDANEK QDGSNFTVCV LPGYTGELCQ

image



Background

Alternative Names Delta and Notch-like epidermal growth factor-related receptor

Note

Product Name: Recombinant Human DNER Protein
Catalog #: PEH3233



For Research Use Only , Not for Diagnostic Use