



Product Name: Recombinant Human RAD50 Protein
Catalog #: PEH3237



Summary

Name	RAD50_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q92878
Host	E.coli
Species	Human
Molecular Mass	70 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LIKAERELSK AEKNSNVETL KMEVISLQNE KADLDRTLRLK LDQEMEQLNH HTTTRTQMEM LTKDKADKDE QIRKIKSRHS DELTSLGYP PNKKQLEDWL HSKSKEINQT RDRLAKLNKE LASSEQNKNH INNELKRKEE QLSSYEDKLF DVCGSQDFES DLDRLKEEIE KSSKQRAMLA GATAVYSQFI TQLTDENQSC CPVCQRVFQT EAELQEVID LQSKLRLAPD KLKSTESSELK KKEKRRDEML GLVPMRQSII DLKEKEIPEL RNKLQNVNRD IQRLKNDIEE QETLLGTIMP EEESAKVCLT DVTIMERFQM ELKDVERKIA QQAALQIGID LDRTVQQVNQ EKQEKQHKLD TVSSKIELNR KLIQDQQEQI QHLKSTTNEL KSEKLQISTN LQRRQQLEEQ TVELSTEVQS LYREIKDAKE QVSPLETTLE KFQQEKEELI NKKNTSNKIA QDKLNDIKEK VKNIHGYMKD IENYIQDGKD DYKKQKETEL NKVIAQLSEC EKHKEKINED MRLMRQDIDT QKIQRWLQD NLTLRKRNEE LKEVEEERKQ HLKEMGQMQV

image



Product Name: Recombinant Human RAD50 Protein
Catalog #: PEH3237



Background

Alternative Names	DNA repair protein RAD50 (hRAD50) (EC 3.6.-.-)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use