


Product Name: Recombinant Human PTPRN2 Protein
Catalog #: PEH3238



Summary

Name	PTPR2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q92932
Host	E.coli
Species	Human
Molecular Mass	11 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	TYSKDLLGQQ PHSEPGAAAF GELQNQMGP SKEEQSLPAG AQEALSDGLQ LEVQPSEEEA RGYIVTDRDP LRPEEGRRLV

image



Background

Alternative Names	Receptor-type tyrosine-protein phosphatase N2 (R-PTP-N2) (EC 3.1.3.-) (EC 3.1.3.48) (Islet cell autoantigen-related protein) (IAR) (ICAAR) (Phogrin) [Cleaved into: IA-2beta60]
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use