


Product Name: Recombinant Human CNDP1 Protein
Catalog #: PEH3246



Summary

Name	CNDP1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	Q96KN2
Host	E.coli
Species	Human
Molecular Mass	63 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	YIDLHQDEFV QTLKEWVAIE SDSVQPVPRF RQELFRMMAV AADTLQRLGA RVASVDMGPQ QLPDGGSLPI PPVILAE LGS DPTKGTVC FY GHLDVQPADR GDGWLTPYV LTEVDGKLYG RGATDNKGPV LAWINAVSAF RALEQDLPVN IKFIIEGMEE AGSVALEELV EKEKDRFFSG VDYIVISDNL WISQRKPAIT YGTRGNSYFM VEVKCRDQDF HSGTFFGILH EPMADLVALL GSLVDSSGHI LVPGIYDEVV PLTEEEINTY KAIHLDLEEY RNSSRVEKFL FDTKEEILMH LWRYP SLIH GIEGAFDEPG TKTVIPGRVI GKFSIRLVPH MNVSAVEKQV TRHLEDVFSK RNSSNKMVVS MTLGLHPWIA NIDDTQYLAA KRAIRTVFGT EPDMIRDGST IPIAKMFQEI VHKS VVLIPL GAVDDGEHSQ

image



Background

Product Name: Recombinant Human CNDP1 Protein
Catalog #: PEH3246



Alternative Names

Beta-Ala-His dipeptidase (EC 3.4.13.20) (CNDP dipeptidase 1) (Carnosine dipeptidase 1) (Glutamate carboxypeptidase-like protein 2) (Serum carnosinase)

Note

For Research Use Only , Not for Diagnostic Use