


Product Name: Recombinant Human ERBIN Protein
Catalog #: PEH3249



Summary

Name	ERBIN_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q96RT1
Host	E.coli
Species	Human
Molecular Mass	67 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	IFTFEKTL EE LYLDANQIEE LPKQLFNCQS LHKLSLPDND LTTLPASIAN LINLRELDVS KNGIQEFPEN IKNCKVLTIV EASVNPISKL PDGFSQLLNL TQLYLND AFL EFLPANFGRL TKLQILELRE NQLKMLPKTM NRLTQLERLD LGSNEFTEVP EVLEQLSGLK EFWMDANRLT FIPGFIGSLK QLT YLDVSKN NIEMVEEGIS TCENLQDLLL SSNSLQQLPE TIGSLKNITT LKIDENQLMY LPDSIGGLIS VEELDCSFNE VEALPSSIGQ LTNLRTFAAD HNYLQQLPPE IGSWKNITVL FLHSNKLETL PEEMGDMQKL KVINLSDNRL KNLPFSFTKL QQLTAMWLS D NQSKPLIPLQ KETDSETQKM VLTNYMFPQQ PRTE DVMFIS DNESFNPSLW EEQRKQRAQV AFECDEDKDE REAPPREGNL KRYPTYP PDE LKNMVKT VQT IVHRLKDEET NEDSGRDLKP HEDQQDINKD VGVKTSESTT TVKSKVDERE KYMIGNSVQK ISEPEAEISP GSLPVTANMK ASEN LKHIVN HDDVFEESEE LSSDEEMKMA

image



Product Name: Recombinant Human ERBIN Protein
Catalog #: PEH3249



Background

Alternative Names Erbin (Densin-180-like protein) (ErbB2-interacting protein) (Protein LAP2)

Note

For Research Use Only , Not for Diagnostic Use