


**Product Name: Recombinant Human PARK7 Protein**  
**Catalog #: PEH3251**



## Summary

|                              |  |
|------------------------------|--|
| <b>Name</b>                  | PARK7_HUMAN  |
| <b>Purity</b>                | Greater than 85% as determined by reducing SDS-PAGE  |
| <b>Endotoxin level</b>       | Please contact with the lab for this information   |
| <b>Tag</b>                   | C-His  |
| <b>Accession #</b>           | Q99497   |
| <b>Host</b>                  | E.coli   |
| <b>Species</b>               | Human  |
| <b>Molecular Mass</b>        | 21 kDa   |
| <b>Buffer</b>                | Phosphate buffered saline  |
| <b>Shipping</b>              | The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.   |
| <b>Stability&amp;Storage</b> | Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles. |
| <b>Sequence</b>              | AKGAEEMETV IPVDVMRRAG IKVTVAGLAG KDPVQCSRDV VICPDASLED<br>AKKEGPYDVV VLPGGNLGAQ NLSESAAVKE ILKEQENRKG LIAAICAGPT<br>ALLAHEIGFG SKVTTHPLAK DKMMNGGHYT YSENVRVEKDG LILTSRGPPT<br>SFEFALAIVE ALNGKEVAAQ           |

## image



## Background

|                          |  |
|--------------------------|--|
| <b>Alternative Names</b> | Parkinson disease protein 7 (Maillard deglycase) (Oncogene DJ1) (Parkinsonism-associated deglycase) (Protein DJ-1) (DJ-1) (Protein/nucleic acid deglycase DJ-1) (EC 3.1.2.-) (EC 3.5.1.-) (EC 3.5.1.124) |
|--------------------------|--|

## Note

**Product Name: Recombinant Human PARK7 Protein**  
**Catalog #: PEH3251**



---

For Research Use Only , Not for Diagnostic Use