


Product Name: Recombinant Human PRRX2 Protein
Catalog #: PEH3253



Summary

Name	PRRX2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q99811
Host	E.coli
Species	Human
Molecular Mass	12 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	CAQARKNFSV SHLLDLEEVA AAGRLAARPG ARAEAREGAA REPSGGSSGS EAAPQDGECP SPGRGSAKR KKKQRRNRTT FNSSQLQALE RVFERTHYPD

image



Background

Alternative Names	Paired mesoderm homeobox protein 2 (Paired-related homeobox protein 2) (PRX-2)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use