


Product Name: Recombinant Human ARHGDIG Protein
Catalog #: PEH3254



Summary

Name	GDIR3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q99819
Host	E.coli
Species	Human
Molecular Mass	23 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EGGPPAVDEV LDEAVPEYRA PGRKSLEIR QLDPDDRSLA KYKRVLLGPL PPAVDPSLPN VQVTRLTLLS EQAPGPVVM D LTGDLAVLKD QVFVLKEGVD YRVKISFKVH REIVSGLKCL HHTYRRGLRV DKTVMVGSY GPSAQEYEFV TPVEEAPRGA LVRGPYLVVS LFTDDDRTHH LSWEWGCLC

image



Background

Alternative Names Rho GDP-dissociation inhibitor 3 (Rho GDI 3) (Rho-GDI gamma)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human ARHGDIG Protein
Catalog #: PEH3254

