


Product Name: Recombinant Human RPP25 Protein
Catalog #: PEH3257



Summary

Name	RPP25_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9BUL9
Host	E.coli
Species	Human
Molecular Mass	21 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MENFRKVRSE EAPAGCGAEG GGPGSGPFAD LAPGAVHMRV KEGSKIRNLM AFATASMAQP ATRAI VFSGC GRATTKTVTC AEILKRRLAG LHQVTRLRYR SVREVWQSLP PGPTQGQTPG EPAASLSVLK NVPGLAILLS KDALDPRQPG YQPPNPHPGP SSPPAAPASK RSLGEPAAAGE GSAKRSQPEP GVADEDQTA

image



Background

Alternative Names	Ribonuclease P protein subunit p25 (RNase P protein subunit p25) (EC 3.1.26.5)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human RPP25 Protein
Catalog #: PEH3257

