


Product Name: Recombinant Human API5 Protein
Catalog #: PEH3259



Summary

Name	API5_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9BZZ5
Host	E.coli
Species	Human
Molecular Mass	32 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DSAINAQLDL CEDEDVSIRR QAIKELPQFA TGENLPRVAD ILTQLLQTDD SAEFNLVNNA LLSIFKMDAK GTLGGLFSQI LQGEDIVRER AIKFLSTKLK TLPDEVLTKE VEELILTESK KVLEDVTGEE FVLFMKILSG LKSLQTVSGR QQLVELVAEQ ADLEQTFNPS DPDCVDRLLQ CTRQAVPLFS KNVHSTRFVT YFCEQVLPNL GTLTPVEGL DIQLEVLKLL AEMSSFCGDM EKLETNLRKL FDKLLEYMPL PPEEAENG

image



Background

Alternative Names	Apoptosis inhibitor 5 (API-5) (Antiapoptosis clone 11 protein) (AAC-11) (Cell migration-inducing gene 8 protein) (Fibroblast growth factor 2-interacting factor) (FIF) (Protein XAGL)
--------------------------	---

Product Name: Recombinant Human API5 Protein
Catalog #: PEH3259



Note

For Research Use Only , Not for Diagnostic Use