


Product Name: Recombinant Human IL22 Protein
Catalog #: PEH3260



Summary

Name	IL22_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	Q9GZX6
Host	E.coli
Species	Human
Molecular Mass	59 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	APISSHCRLD KSNFQQPYIT NRTFMLAKEA SLADNNTDVR LIGEKLFHGV SMSERCYLMK QVLNFTLEEV LFPQSDRFQP YMQEVVPFLA RLSNRLSTCH IEGDDLHIQR NVQKLKDTVK KLGESGEIKA IGELDLLFMS LRNACI

image



Background

Alternative Names	Interleukin-22 (IL-22) (Cytokine Zcyto18) (IL-10-related T-cell-derived-inducible factor) (IL-TIF)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use