


Product Name: Recombinant Human AZI2 Protein
Catalog #: PEH3264



Summary

Name	AZI2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9H6S1
Host	E.coli
Species	Human
Molecular Mass	23 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	TSSVGREQVN KAYHAYREVC IDRDNLKSKL DKMNKDNSES LKVLNEQLQS KEVELLQLRT EVETQQVMRN LNPPSSNWEV EKLSCDLKIH GLEQELELMR KECSDLKIEL QKAKQTDPYQ EDNLKSRDLQ KLSISSDNMQ HAYWELKREM SNLHLVTQVQ AELLRKLKTS TAIKKACAPV GCSEDLGRDS

image



Background

Alternative Names 5-azacytidine-induced protein 2 (NF-kappa-B-activating kinase-associated protein 1) (Nak-associated protein 1) (TILP)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human AZI2 Protein
Catalog #: PEH3264

