


Summary

Name	MUTYH_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9UIF7
Host	E.coli
Species	Human
Molecular Mass	31 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	RPGDFNQAAM ELGATVCTPQ RPLCSQCPVE SLCRARQRVE QEQLLAGSLS SGSPDVEECA PNTGQCHLCL PPSEPWDQTL GVVNFPRKAS RKPPRESSA TCVLEQPGAL GAQILLVQRP NSGLLAGLWE FPSVTWEPSE QLQRKALLQE LQRWAGPLPA THLRHLGEVV HTFSHIKLTQ QVYGLALEGQ TPVTTVPPGA RWLTQEEFHT AAVSTAMKKV FRVYQGQQPG TCMGSKRSQV SSPCSRKKPR MGQQVLDNFF

image



Background

Alternative Names Adenine DNA glycosylase (EC 3.2.2.31) (MutY homolog) (hMYH)

Product Name: Recombinant Human MUTYH Protein
Catalog #: PEH3279



Note

For Research Use Only , Not for Diagnostic Use