


Product Name: Recombinant Human MBD1 Protein
Catalog #: PEH3281



Summary

Name	MBD1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9UIS9
Host	E.coli
Species	Human
Molecular Mass	11 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	VASKKRKKPS RPAKTRKRQV GPQSGEVRKE APRDETKADT DTAPASFPAP GCCENCGISF SGDGTQRQRL KTLCKDCRAQ

image



Background

Alternative Names	Methyl-CpG-binding domain protein 1 (CXXC-type zinc finger protein 3) (Methyl-CpG-binding protein MBD1) (Protein containing methyl-CpG-binding domain 1)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use