


Product Name: Recombinant Human GGA1 Protein
Catalog #: PEH3282



Summary

Name	GGA1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9UJY5
Host	E.coli
Species	Human
Molecular Mass	16 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LLKSSHPEDL RAANKLIKEM VQEDQKRMEK ISKRVNAIEE VNNNVKLLTE MVMSHSQGGA AAGSSEDLMK ELYQRCERMR PTLFRLASDT EDNDEALAEI LQANDNLTQV INLYKQLVRG

image



SDS



WB

Background

Alternative Names	ADP-ribosylation factor-binding protein GGA1 (Gamma-adaptin-related protein 1) (Golgi-localized, gamma ear-containing, ARF-binding protein 1)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use