


Product Name: Recombinant Human COL17A1 Protein
Catalog #: PEH3286



Summary

Name	COHA1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9UMD9
Host	E.coli
Species	Human
Molecular Mass	25 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GTEVTERIVT ETVTTRLTSL PPKGTSNGY AKTASLGGS RLEKQSLTHG SSGYINSTGS TRGHASTSSY RRAHSPASTL PNSPGSTFER KTHVTRHAYE GSSSGNSSPE YPRKEFASSS TRGRSQTRES EIRVRLQSAS PSTRWTELDD VKRLLKGSRS ASVSPTRNSS NTLPIPKKGT VETKIVTASS QSVSGTYDAT ILDANLPSHV WSSTLPAGSS

image



Background

Alternative Names	Collagen alpha-1(XVII) chain (180 kDa bullous pemphigoid antigen 2) (Bullous pemphigoid antigen 2) [Cleaved into: 120 kDa linear IgA disease antigen (120 kDa linear IgA dermatosis antigen) (Linear IgA disease antigen 1) (LAD-1); 97 kDa linear IgA disease antigen (97 kDa linear IgA bullous dermatosis antigen)]
--------------------------	--

Product Name: Recombinant Human COL17A1 Protein
Catalog #: PEH3286



(97 kDa LAD antigen) (97-LAD) (Linear IgA bullous disease antigen of 97 kDa)
(LABD97)]

Note

For Research Use Only , Not for Diagnostic Use