


Product Name: Recombinant Human SIGLEC9 Protein
Catalog #: PEH3294



Summary

Name	SIGL9_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9Y336
Host	E.coli
Species	Human
Molecular Mass	11 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	AAGVGDTGIE DANA VRGSAS QGPLTEPWAE DSPPDQPPPA SARSSVGEGE LQYASLSFQM VKPWDSRGQE ATDTEYSEIK

image



Background

Alternative Names	Sialic acid-binding Ig-like lectin 9 (Siglec-9) (CDw329) (Protein FOAP-9) (CD antigen CD329)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use