


Product Name: Recombinant Human HYOU1 Protein
Catalog #: PEH3298



Summary

Name	HYOU1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9Y4L1
Host	E.coli
Species	Human
Molecular Mass	26 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SEAQKPSEKA EAGPEGVAPA PEGEKKQKPA RKRRMVEEIG VELVVLDLPD LPEDKLAQSV QKLQDLTLRD LEQQEREKAA NSLEAFIFET QDKLYQPEYQ EVSTEEQREE ISGKLSAAST WLEDEGVGAT TVMLKEKLA E LRLCQGLFF RVEERKKWPE RLSALDNLLN HSSMFLKGAR LIPEMDQIFT EVEMTTLEKV INETWAWKNA

image



Background

Alternative Names	Hypoxia up-regulated protein 1 (150 kDa oxygen-regulated protein) (ORP-150) (170 kDa glucose-regulated protein) (GRP-170)
--------------------------	---

Note

Product Name: Recombinant Human HYOU1 Protein
Catalog #: PEH3298



For Research Use Only , Not for Diagnostic Use