


Product Name: Recombinant Human ETV7 Protein
Catalog #: PEH3300



Summary

Name	ETV7_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9Y603
Host	E.coli
Species	Human
Molecular Mass	13 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ADCRLLDYV YQLLLDTRYE PYIKWEDKDA KIFRVVDPNG LARLWGNHKN RVNMTYEKMS RALRHYYKLN IIKKEPGQKL LFRFLKTPGK

image



Background

Alternative Names Transcription factor ETV7 (ETS translocation variant 7) (ETS-related protein Tel2) (Tel-related Ets factor) (Transcription factor Tel-2)

Note

For Research Use Only , Not for Diagnostic Use