


**Product Name: Recombinant Human GFAP Protein**  
**Catalog #: PEH3304**



## Summary

<b>Name</b>	GFAP_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-His
<b>Accession #</b>	P14136
<b>Host</b>	E.coli
<b>Species</b>	Human
<b>Molecular Mass</b>	39 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	GFKETRASER AEMMELNDRF ASYIEKVRFL EQQNKALAAE LNQLRAKEPT KLADVYQAEI RELRLRLDQL TANSARLEVE RDNLAQDLAT VRQKLQDET LRLEAENNLA AYRQEADEAT LARLDLERKI ESLEEEIRFL RKIHEEEVRE LQEQLARQQV HVELDVAKPD LTAALKEIRT QYEAMASSNM HEAEWYRSK FADLTDAAR NAELLRQAKH EANDYRRQLQ SLTCDLESLR GTNESLERQM REQEERHVRE AASYQEALAR LEEEGQSLKD EMARHLQEYQ DLLNVKLALD IEIATYRKLL

## image



## Background

**Alternative Names** Glial fibrillary acidic protein (GFAP)

**Product Name: Recombinant Human GFAP Protein**  
**Catalog #: PEH3304**

---



**Note**

For Research Use Only , Not for Diagnostic Use