


Product Name: Recombinant Human IFNA1 Protein
Catalog #: PEH3310



Summary

Name	IFNA1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P01562
Host	E.coli
Species	Human
Molecular Mass	20 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	CDLPETHSLD NRRTLMLLAQ MSRISPSSCL MDRHDFGFQK EEFDGNQFQK APAISVLHEL IQQIFNLFTT KDSSAAWDED LLDKFCTELY QQLNDLEACV MQEERVGETP LMNADSILAV KKYFRRITLY LTEKKYSPCA WEVVRAEIMR SLSLSTNLQE RLRRKE

image



Background

Alternative Names Interferon alpha-1/13 (IFN-alpha-1/13) (Interferon alpha-D) (LeIF D)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human IFNA1 Protein
Catalog #: PEH3310

