


Product Name: Recombinant Human GPR15LG Protein
Catalog #: PEH3322



Summary

Name	GP15L_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-GST
Accession #	Q6UWK7
Host	E.coli
Species	Human
Molecular Mass	32 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KRRPAKAWSG RRRLCCHR V PSPNSTNLKG HHVRLCKPCK LEPEPRLWVV PGALPQV

image



Background

Alternative Names	Protein GPR15L (Antimicrobial peptide with 57 amino acid residues) (AP-57) (Antimicrobial peptide-57) (Colon-derived SUSD2 binding factor) (CSBF) (Secreted protein C10orf99)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use