


**Product Name: Recombinant Human RHNO1 Protein**  
**Catalog #: PEH3323**



## Summary

|                              |                                                                                                                                                                                                                 |
|------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Name</b>                  | RHNO1_HUMAN                                                                                                                                                                                                     |
| <b>Purity</b>                | Greater than 85% as determined by reducing SDS-PAGE                                                                                                                                                             |
| <b>Endotoxin level</b>       | Please contact with the lab for this information                                                                                                                                                                |
| <b>Tag</b>                   | C-His                                                                                                                                                                                                           |
| <b>Accession #</b>           | Q9BSD3                                                                                                                                                                                                          |
| <b>Host</b>                  | E.coli                                                                                                                                                                                                          |
| <b>Species</b>               | Human                                                                                                                                                                                                           |
| <b>Molecular Mass</b>        | 16 kDa                                                                                                                                                                                                          |
| <b>Buffer</b>                | Lyophilized from Phosphate buffered saline                                                                                                                                                                      |
| <b>Shipping</b>              | The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.                                                                                              |
| <b>Stability&amp;Storage</b> | Store at $\leq -70^{\circ}\text{C}$ , stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles. |
| <b>Sequence</b>              | HSSRKPTTSK FPHLTFESQ SSSSETLGIP LIRECPSESE KDVSRRPLVP<br>VLSPQSCGNM SVQALQSLPY VFIPPDQTP ESSSVKEELI PQDQKENSLL<br>SCTLHTGTPN SPEPGPVLVK                                                                         |

## image



## Background

|                          |                                                                                                             |
|--------------------------|-------------------------------------------------------------------------------------------------------------|
| <b>Alternative Names</b> | RAD9, HUS1, RAD1-interacting nuclear orphan protein 1 (RAD9, RAD1, HUS1-interacting nuclear orphan protein) |
|--------------------------|-------------------------------------------------------------------------------------------------------------|

## Note

For Research Use Only , Not for Diagnostic Use