


Product Name: Recombinant Human ATG9A Protein
Catalog #: PEH3325



Summary

Name	ATG9A_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	Q7Z3C6
Host	E.coli
Species	Human
Molecular Mass	26 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	AQAEPERHVW HRRESDESGE SAPDEGGEGA RAPQSIPRSA SYPCAAPRPG APETTALHGG FQRRYGGITD PGTVPRVPSH FSRLPLGGWA EDGQSASRHP EPVPEEGSED

image



Background

Alternative Names Autophagy-related protein 9A (APG9-like 1) (mATG9)

Note

For Research Use Only , Not for Diagnostic Use