


**Product Name: Recombinant Human TRPM3 Protein**  
**Catalog #: PEH3327**



## Summary

<b>Name</b>	TRPM3_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-His
<b>Accession #</b>	Q9HCF6
<b>Host</b>	E.coli
<b>Species</b>	Human
<b>Molecular Mass</b>	20 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	DYGLKLFITD DELKKVHDFE EQCIEEYFRE KDDRFNSSND ERIRVTSERV ENMSMRLEEV NEREHSMKAS LQTVDIRLAQ LEDLIGRMAT ALERLTGLER AESNKIRSRT SSDCTDAAYI VRQSSFNSQE GNTFKLQESI DPAGEETMSP TSPTLMPRMR

## image



## Background

<b>Alternative Names</b>	Transient receptor potential cation channel subfamily M member 3 (Long transient receptor potential channel 3) (LTrpC-3) (LTrpC3) (Melastatin-2) (MLSN2)
--------------------------	--

## Note

**Product Name: Recombinant Human TRPM3 Protein**  
**Catalog #: PEH3327**



---

For Research Use Only , Not for Diagnostic Use