


Product Name: Recombinant Human TRPM7 Protein
Catalog #: PEH3329



Summary

Name	TRPM7_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q96QT4
Host	E.coli
Species	Human
Molecular Mass	17 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DQETVCSKAT EGDNTEFGAF VGHRDSMDLQ RFKETSNIKILSNNNTSEN TLKRVSSLAG FTDCHRTSIP VHSKQAEKIS RRPSTEDTHE VDSKAALIPD WLQDRPSNRE MPSEEGTLNG LTSPFKPAMD

image



Background

Alternative Names	Transient receptor potential cation channel subfamily M member 7 (EC 2.7.11.1) (Channel-kinase 1) (Long transient receptor potential channel 7) (LTrpC-7) (LTrpC7)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human TRPM7 Protein
Catalog #: PEH3329

