


Product Name: Recombinant Human TP63 Protein
Catalog #: PEH3337



Summary

Name	P63_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9H3D4
Host	E.coli
Species	Human
Molecular Mass	12 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DRKAEDSIR KQQVSDSTKN GDGTKRPFRQ NTHGIQMTSI KKRRSPDDEL LYLPVRGRET YEMLLKIKES LELMQYLPQH TIETYRQQQQ

image



Background

Alternative Names Tumor protein 63 (p63) (Chronic ulcerative stomatitis protein) (CUSP) (Keratinocyte transcription factor KET) (Transformation-related protein 63) (TP63) (Tumor protein p73-like) (p73L) (p40) (p51)

Note

For Research Use Only , Not for Diagnostic Use