


Product Name: Recombinant Human SNRNP70 Protein
Catalog #: PEH3342



Summary

Name	RU17_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P08621
Host	E.coli
Species	Human
Molecular Mass	23 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ALFAPRDIPI YLPPLEKLPH EKHHNQPYCG IAPYIREFED PRDAPPPTRA ETREERMERK RREKIERRQQ EVETELKMWD PHNDPNAQGD AFKTLFVARV NYDTTESKLR REFEVYGPIK RIHMOVYSKRS GKPRGYAFIE YEHERDMHSA YKHADGKKID GRRVLVDVER

image



Background

Alternative Names U1 small nuclear ribonucleoprotein 70 kDa (U1 snRNP 70 kDa) (U1-70K) (snRNP70)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human SNRNP70 Protein
Catalog #: PEH3342

