


Product Name: Recombinant Human PTPRN Protein
Catalog #: PEH3343



Summary

Name	PTPRN_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q16849
Host	E.coli
Species	Human
Molecular Mass	22 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	VGADIKKTME GPVEGRDTAE LPARTSPMPG HPTASPTSSE VQQVPSPVSS EPPKAARPPV TPVLEKKSPP LGQSQPTVAG QPSARPAEE YGYIVTDQKP LSLAAGVKLL EILAEHVHMS SGSEFINISVV GPALTFRIRH NEQNLSLADV TQQAGLVKSE LEAQTGLQLL QTGVGQREEA AAVLPQTAHS

image



Background

Alternative Names	Receptor-type tyrosine-protein phosphatase-like N (R-PTP-N) (Islet cell antigen 512) (ICA 512) (Islet cell autoantigen 3) (PTP IA-2) [Cleaved into: ICA512-N-terminal fragment (ICA512-NTF); ICA512-transmembrane fragment (ICA512-TMF); ICA512-cleaved cytosolic fragment (ICA512-CCF)]
--------------------------	--

Product Name: Recombinant Human PTPRN Protein
Catalog #: PEH3343



Note

For Research Use Only , Not for Diagnostic Use