


Product Name: Recombinant Human PTPRN Protein
Catalog #: PEH3344



Summary

Name	PTPRN_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q16849
Host	E.coli
Species	Human
Molecular Mass	35 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GIFGDHPGHS YGDLPGPSA QLFQDSGLLY LAQELPAPSR ARVPRLPEQG SSSRAEDSPE GYEKEGLGDR GEKPASPAVQ PDAALQRLAA VLAGYGVELR QLTPEQLSTL LTLQLLPKG AGRNPGGVVN VGADIKKTME GPVEGRDTAE LPARTSPMPG HPTASPTSSE VQQVPSPVSS EPPKAARPPV TPVLEKKSP LGQSQPTVAG QPSARPAEE YGYIVTDQKP LSLAAGVKLL EILAEHVHMS SGSFINISVV GPALTFRIRH NEQNLSLADV TQQAGLVKSE LEAQTGLQL QTGVGQREEA

image



Background

Alternative Names	Receptor-type tyrosine-protein phosphatase-like N (R-PTP-N) (Islet cell antigen 512) (ICA 512) (Islet cell autoantigen 3) (PTP IA-2) [Cleaved into:
--------------------------	---

Product Name: Recombinant Human PTPRN Protein
Catalog #: PEH3344



ICA512-N-terminal fragment (ICA512-NTF); ICA512-transmembrane
fragment (ICA512-TMF); ICA512-cleaved cytosolic fragment (ICA512-CCF)]

Note

For Research Use Only , Not for Diagnostic Use