


Product Name: Recombinant Human NCAM1 Protein
Catalog #: PEH3349



Summary

Name	NCAM1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	P13591
Host	E.coli
Species	Human
Molecular Mass	34 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QVDIVPSQGE ISVGESKFFL CQVAGDAKDK DISWFSPNGE KLTPNQQRIS VWVNDDSSST LTIYNANIDD AGIYKCVVTG EDGSESEATV NVKIFQKLMF KNAPTPQEFR EGEDAVIVCD VVSSLPTII WKHKGRDVIL KKDVRFIVLS NNYLQIRGIK KTDEGTYRCE GRILARGEIN

image



Background

Alternative Names Neural cell adhesion molecule 1 (N-CAM-1) (NCAM-1) (CD antigen CD56)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human NCAM1 Protein
Catalog #: PEH3349

