


Product Name: Recombinant Human CD22 Protein
Catalog #: PEH3351



Summary

Name	CD22_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P20273
Host	E.coli
Species	Human
Molecular Mass	13 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LQENSSGQSF FVRNKKVRRR PLSEGPLSLG CYNPMMEDGI SYTTLRFPPEM NIPRTGDAES SEMQRPPPDC DDTVTYSALH KRQVGDYENV

image



Background

Alternative Names B-cell receptor CD22 (B-lymphocyte cell adhesion molecule) (BL-CAM) (Sialic acid-binding Ig-like lectin 2) (Siglec-2) (T-cell surface antigen Leu-14) (CD antigen CD22)

Note

For Research Use Only , Not for Diagnostic Use