


Summary

Name	KIT_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	P10721
Host	E.coli
Species	Human
Molecular Mass	33 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PGESPSPSIH PGKSDLIVRV GDEIRLLCTD PGFVKWTFEI LDETENENKQN EWITEKAEAT NTGKYTCTNK HGLSNSIYVF VRDPAKFLV DRSLYGKEDN DTLVRCPDLD PEVTNYSLKG CQGKPLPKDL RFIPDPKAGI MIKSVKRAYH RLCLHCSVDQ EGKSVLSEKF

image



Background

Alternative Names	Mast/stem cell growth factor receptor Kit (SCFR) (EC 2.7.10.1) (Piebald trait protein) (PBT) (Proto-oncogene c-Kit) (Tyrosine-protein kinase Kit) (p145 c-kit) (v-kit Hardy-Zuckerman 4 feline sarcoma viral oncogene homolog) (CD antigen CD117)
--------------------------	---

Product Name: Recombinant Human KIT Protein
Catalog #: PEH3356



Note

For Research Use Only , Not for Diagnostic Use