



Product Name: Recombinant Human EZH2 Protein
Catalog #: PEH3366



Summary

Name	EZH2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q15910
Host	E.coli
Species	Human
Molecular Mass	20 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GPQCYQHLEG AKEFAAALTA ERIKTPPKRP GGRRRGRLPN NSSRPSTPTI NVLESKDTDS DREAGTETGG ENNDKEEEK KDETSSSEA NSRCQTPIKM KPNIEPPENV EWSGAEASMF RVLIGTYIDN FCAIARLIGT KTCRQVYEFR VKESSIAPA

image



Background

Alternative Names Histone-lysine N-methyltransferase EZH2 (EC 2.1.1.356) (ENX-1) (Enhancer of zeste homolog 2) (Lysine N-methyltransferase 6)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human EZH2 Protein
Catalog #: PEH3366

