


Product Name: Recombinant Human ITGB1 Protein
Catalog #: PEH3371



Summary

Name	ITB1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P05556
Host	E.coli
Species	Human
Molecular Mass	13 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	NAKSCGECIQ AGPNCGWCTN STFLQEGMPT SARCDLEAL KKKGCPPDDI ENPRGSKDIK KNKNVTNRSK GTAEKLPED ITQIQPQQLV LRLRSGEPQT

image



Background

Alternative Names Integrin beta-1 (Fibronectin receptor subunit beta) (Glycoprotein IIa) (GPIIA) (VLA-4 subunit beta) (CD antigen CD29)

Note

For Research Use Only , Not for Diagnostic Use