


Product Name: Recombinant Human CD44 Protein
Catalog #: PEH3375



Summary

Name	CD44_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P16070
Host	E.coli
Species	Human
Molecular Mass	20 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GYRQTPKEDS HSTTGTAAS AHTSHPMQGR TTPSPEDSSW TDFNPISHP MGRGHQAGRR MDMDSSHSIT LQPTANPNTG LVEDLDRTGP LSMTTQQSNS QSFSTSHEGL EEDKDHPPTS TLSSNRNDV TGGRRDPNHS EGSTTLLEGY TSHYPHTKES

image



Background

Alternative Names	CD44 antigen (CDw44) (Epican) (Extracellular matrix receptor III) (ECMR-III) (GP90 lymphocyte homing/adhesion receptor) (HUTCH-I) (Heparan sulfate proteoglycan) (Hermes antigen) (Hyaluronate receptor) (Phagocytic glycoprotein 1) (PGP-1) (Phagocytic glycoprotein I) (PGP-I) (CD antigen CD44)
--------------------------	--

Product Name: Recombinant Human CD44 Protein
Catalog #: PEH3375



Note

For Research Use Only , Not for Diagnostic Use