


Summary

Name	VEGFA_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	P15692
Host	E.coli
Species	Human
Molecular Mass	58 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MRCGGCCNDE GLECVPTES NITMQIMRIK PHQGQHIGEM SFLQHNKCEC RPKKDRARQE KKSVRGKGGK QKRKRKKSRY KSWSVYVGAR CCLMPWSLPG PHPCGPCSER RKHLFVQDPQ TCKCCKNTD

image



Background

Alternative Names	Vascular endothelial growth factor A (VEGF-A) (Vascular permeability factor) (VPF)
--------------------------	------------------------------------------------------------------------------------

Note

For Research Use Only , Not for Diagnostic Use