


Product Name: Recombinant Human FOXP3 Protein
Catalog #: PEH3386



Summary

Name	FOXP3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9BZS1
Host	E.coli
Species	Human
Molecular Mass	12 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	FKFHNMRPPF TYATLIRWAI LEAPEKQRTL NEIYHWFTRM FAFFRNHPAT WKNAIRHNLS LHKCFVRVES EKGAVWTVDE

image



Background

Alternative Names	Forkhead box protein P3 (Scurfin) [Cleaved into: Forkhead box protein P3, C-terminally processed; Forkhead box protein P3 41 kDa form]
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use